



proven high-performance automated liquid handling

apricot ipp

features

- Proven high-performance automated liquid handling
- Control via easy to use interface – Apple iPod® or Apple iPad mini® controller
- EZ-Load® Tip Technology – Load tips quickly and easily, without tools!
- Automated plate loading via motorized elevator
- 2 or 3-position shuttle /programmable serial dilution
- Small footprint maximizes bench space and fits suitable hoods

applications

- Plate replication
- Plate reformatting
- Reagent addition
- Compound addition
- Serial dilution by column

models

- 96L 0.5–125µL
- 96H 1–500µL
- 96E 5–1000µL
- 384 0.5–125µL

Ask us about custom configurations & accessories.

specifications	• Dimensions:	711mm x 483mm x 607mm (9"W x 15"D x 14"H)
	• Weight:	45-52 lbs. (45 kg)
	• Power:	100-240 VAC 50/60Hz
	• Operating Temp	50-100° F (10-37° C)